
Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that the our products

e-motion M10 and e-motion M12

including the nickel metal hydridebatteries(Special Provision A123) fulfil the requirements of

IATA chapter 2.3.2.2

and therefore may be carried installed onthe wheelchair as checked baggage without restriction.

Condition precedent to meet the above listed requirement is the transportation in a way that an unintentionally switching on and damage by other luggage is avoided.**The batteries must be removed for transportation** (protected from damage and from short circuit) and must be carried in in strong, rigid packaging in thecargo compartment.

It is allowed to carry spare batteries in the passenger cabin according to IATA chapter 2.3.A.

The pilot-in-command must be informed of the location of thebattery.

We strongly recommend obtaining the approval by corresponding airlinecompany before commencement date respectively departure date.

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur onan annual basis this certificate is valid until **31st December 2014**.

Albstadt, 02th of January 2014

Ulrich Alber GmbH



Ralf Ledda
Managing Director